Basics of SPECT, PET and PET/CT Imaging

Talk presented at the:

Joint Advanced Student School in the State University of St. Petersburg

April 2nd-12th, 2006

Marisa Bazañez-Borgert











Outline

Radiation protection principles

Nuclear medicine

Tracer principle

Radionuclid production

Detection

Gamma camera

Tomographic techniques: SPECT/PET

CT

PET/CT



Some useful units

Activity unit: **1 Bequerel = 1 Decay/s**

Dosis unit:

1 Gray = 1 J/kg 1 Gy warm 1I Water 0.00025 °C

Siever (Sv) refers to equivalent dose in which we measure not only the deposited energy of the particle per unit mass but the effect of this radiation in the tissue. For the photons and energies used in diagnostic imaging:

1 Gray = 1 Sv





Effective dose for different imaging techniques

Image method	Effective Dose (mSv)
Pelvis or hips (x-ray)	0,83 mSv
СТ	10 mSv
PET (FDG)	7 mSv
PET/CT	14 mSv





Effects of whole body exposition

0 -0,25 Sv	no acute effects
1-2 Sv	nausea, vomitus
2-3 Sv	nausea, vomitus on the 1st day recovery after 3 months
3-6 Sv	nausea, vomitus after some hours, hair lose, internal haemorrhages after 2 weeks
6-8 Sv	nausea, vomitus and diarrhea after some hours, mental confusion
50-100 Sv	acute toxic effects, cramps, internal haemorrhages
> 1000 Sv	nerve system damage, dead within seconds



Radiation's effects

Physical: Ionisation and excitation

Physical-chemical: Radicals

Chemical: Reaction of the radicals with bio-molecules

Biological: Cellular dead, mutations





Deterministic effects

They have a threshold dose and their severity depends on the received dose

Stochastic effects

There is a probability of occurrence depending on the dose











Radiopharmaceuticals

are compounds of biomedical interest which have attached radionuclides as labels that allow us to quantify specific physiological processes.

Selection

of pharmaceutical based on organ-specific question

Labelling

of pharmaceutical with radioactive isotopes

Radiopharmaceuticals should not disturb the process under investigation





Non invasive determination of physiologic processes

Tracer principle:

Radiopharmaceuticals

are distributed, metabolized, and excreted according to their chemical structure



Detect and display as:

- Images
- Numerical data
- Time-activity curves





Nuclear Stability













			Pro	duct	Deca	ay Mode		oundance o otope (%)	F			
	Particle acceleration: Cyclotron		¹¹ C ¹³ N ¹⁵ O ¹⁸ F		β	+, EC	9	9.6				
						β+	9	9.8				
					β+ β+, EC		99.6 0.20		ΔΦ			
			⁶⁷ Ga		(EC,γ)		18.8		-			
			¹¹¹ In ¹²³ I ²⁰¹ TI		(EC,γ) (EC,γ) (EC,γ)		48.2 2.6 13.2		Time	Time		
	Radionuclide	nuclide Parent		Halt	flife	Decay mode	Daugther	Half life	Decay mode	Product		
	generators	⁹⁹ Mo		66,0)2 h	β-	^{99m} Tc	6,02 h	Isomer	⁹⁹ Tc		
	68G	¹⁸⁸ W	69,4		ŀd	β-	¹⁸⁸ Re	16,9 h	β-	¹⁸⁸ Os		
		⁶⁸ Ge		288	d	EC	⁶⁸ Ga	68,3 min	β+, EC	⁶⁸ Zn		
		⁸² Sr		25 c	ł	EC	⁸² Rb	1,25 min	β+, EC	⁸² Kr		





Interaction of Radiation with Matter







Interaction of Radiation with Matter



Part of the γ -Energy is passed to the electron the γ -Quantum changes its Energy and direction





Probability of Interaction

Probability of interaction versus photon energy (for absorbers of different atomic numbers)





Interaction

May pass through May interact May interact no interaction lose energy & change direction transfer all its energy & disappear



$$\mu_m(mix) = \mu_{m1}f_1 + \mu_{m2}f_2 + \dots$$

 $\mu_m = \tau + \sigma + \kappa$

- $_{.\tau}$ is the part of μ_m due to photoelectric effect
- σ is the part of μ_m due to Compton Effect
- κ is the part of μ_m due to pair production
- μ_m Absorption coefficient involves both absorption and scattering processes.







Ionised electron is Electron-Hole Pair Activator de-excites by emission of Photons (10⁻⁷ sec)





Scintillation Crystals

	Nal	BGO	LSO	GSO
	Nal:TI	Bi4Ge3O12	Lu2SiO5:Ce	Gd2SiO5:Ce
Density (g/cm ³)	3,67	7,1	7,4	6,7
Effective atomic No.	51	73	65	59
photons/MeV	41000	9000	26000	8000
Decay time (ns)	230	300	40	60
Peak emission (nm)	415	480	420	430
Index of refraction	1.85	2.15	1.82	1.85





Scintillation Crystals: Ideal Properties

- High photoconversion
- High light output
- Good stopping power
- Fast response
- High optical light transmission and no selfabsorption
- Index of refraction close to that of the PMT
- Nonhygroscopic and robust



Scintillators







Photomultiplier Tubes







Scintillation Crystal







The Gamma Camera







Principle of the Gamma Camera









Energy Spectrum







Collimators

- Diverging Hole
 - Large FOV; minifies large objects
 - Poor sensitivity, resolution
- Parallel Hole
 - No magnification
 - FOV, sensitivity constant with distance
 - Resolution limited by camera intrinsic resolution
- Pinhole, Converging Hole
 - Magnification of small objects
 - Sensitivity and resolution good at close distances





The Objective of the Gamma Camera is to get an image of the distribution of the activity









Behind the scene...





SPECT









Tomographic Images







Measured data grouped into projections: SINOGRAM **90°** Projection at 90° . **45°** ° Projection at 0° -45° S





Image Reconstruction Methods

Simple Backprojection Filtered Backprojection Fourier Transformation

Iterative Algorithms













Electronic Collimation






PET Acquisition Modes







PET Coincidence Types





True Coincidences

Scattered within FOV

Scattered outside FOV

Random Coincidences

High Crystal Efficiency Good Energy Resolution Short Coincidence Window











PET Clinical Applications







PET Clinical Applications

Definition of tumor extension



tumor volume: 4 ml









PET Clinical Applications Glucose Metabolism in the Brain Surface Projections











CT







CT





















CT–based attenuation coefficients







Research at NMK in KRI



Respiratory Movement









... and the winner is... Thanks to our Russian colleagues for the great time together





Thank you very much for your attention

Hey, there is some time left today and It's Tuesday ...

let's go BOWLING!!!!!



